



EPSA[®] *safety first*

**CATHETER FOR
URETEROCUTANEOSTOMY**

CATHETER FOR URETEROCUTANEOSTOMY

CODE: 9401 (26 cm) – 9402 (32 cm)

1 PUR catheter, flexible, multi-fenestrated, with antiperistaltic anchorage fins, and disc for skin fixation.

STERILE, SINGLE PACK, SINGLE USE.



Available Sizes

3,5 / 4,2 mm

Length: 26 cm. / 32 cm.

Recommendation for use

Recommendation of a professional on the use of the technology.

The Catheter for Ureterocutaneostomy is a therapeutic device designed to calibrate the terminal sector of a ureterocutaneostomy.

(Ureterocutaneostomy: surgical procedure consisting of a communication of the ureter with the outside through a cutaneous hole).

It is indicated in those patients in whom an advanced malignant pathological process (bladder, prostate or rectal tumors with infiltrative growth) does not allow the restoration of urinary transit, or a short life expectancy does not justify the performance of internal derivative procedures (creation of reservoir or neo-bladder, etc.).

